

AGXT2L2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16978	SYNONYMS
SIZE	100 ug	AGXT2L2; PP9286; 5-phosphohydroxy-L-lysine
HOST	Rabbit	phospho-lyase; Alanine--glyoxylate aminotransferase 2-like 2

PRODUCT INFORMATION

GENE NAME	AGXT2L2
PROTEIN NAME	5-phosphohydroxy-L-lysine phospho-lyase
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human AGXT2L2. at AA range: 310-390

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.