

Mnk1 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17009	SYNONYMS
SIZE	100 ug	MKNK1; MNK1; MAP kinase-interacting
HOST	Rabbit	serine/threonine-protein kinase 1; MAP kinase signal-integrating kinase 1; MAPK signal-integrating kinase 1; Mnk1

PRODUCT INFORMATION

GENE NAME	MKNK1
PROTEIN NAME	MAP kinase-interacting serine/threonine-protein kinase 1
IMMUNOGEN	Synthesized peptide derived from human Mnk1 around the non-phosphorylation site of T255.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.