

Olfactory receptor 6A2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17011	SYNONYMS
SIZE	100 ug	
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	OR6A2
PROTEIN NAME	Olfactory receptor 6A2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Olfactory receptor 6A2. at AA range: 30-110

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.