

PEA-15 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17038	SYNONYMS
SIZE	100 ug	PEA15; Astrocytic phosphoprotein PEA-15; 15 kDa
HOST	Rabbit	phosphoprotein enriched in astrocytes; Phosphoprotein enriched in diabetes; PED

PRODUCT INFORMATION

GENE NAME	PEA15
PROTEIN NAME	Astrocytic phosphoprotein PEA-15
IMMUNOGEN	Synthesized peptide derived from human PEA-15 around the non-phosphorylation site of S116.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat,Monkey
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.