

Ksr-1 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17071	SYNONYMS
SIZE	100 ug	KSR1; KSR; Kinase suppressor of Ras 1
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	KSR1
PROTEIN NAME	Kinase suppressor of Ras 1
IMMUNOGEN	Synthesized peptide derived from human Ksr-1 around the non-phosphorylation site of S392.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse
TESTED APPLICATIONS	WB, IHC, IF, ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.