

eIF2 β Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17087	SYNONYMS
SIZE	100 ug	EIF2S2; EIF2B; Eukaryotic translation initiation factor 2
HOST	Rabbit	subunit 2; Eukaryotic translation initiation factor 2 subunit beta; eIF-2-beta

PRODUCT INFORMATION

GENE NAME	EIF2S2
PROTEIN NAME	Eukaryotic translation initiation factor 2 subunit 2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human eIF2 β . at AA range: 30-110

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.