

EDG-4 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17158	SYNONYMS
SIZE	100 ug	LPAR2; EDG4; LPA2; Lysophosphatidic acid receptor
HOST	Rabbit	2; LPA receptor 2; LPA-2; Lysophosphatidic acid receptor Edg-4

PRODUCT INFORMATION

GENE NAME	LPAR2
PROTEIN NAME	Lysophosphatidic acid receptor 2
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human EDG-4. at AA range: 250-330

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Monkey
TESTED APPLICATIONS	WB,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.