

Ribosomal Protein L23 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17207	SYNONYMS
SIZE	100 ug	RPL23; 60S ribosomal protein L23; 60S ribosomal
HOST	Rabbit	protein L17

PRODUCT INFORMATION

GENE NAME	RPL23
PROTEIN NAME	60S ribosomal protein L23
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Ribosomal Protein L23. at AA range: 30-110

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat,Monkey
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.