

Gab 2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17292	SYNONYMS
SIZE	100 ug	GAB2; KIAA0571; GRB2-associated-binding protein 2;
HOST	Rabbit	GRB2-associated binder 2; Growth factor receptor bound protein 2-associated protein 2; pp100

PRODUCT INFORMATION

GENE NAME	GAB2
PROTEIN NAME	GRB2-associated-binding protein 2
IMMUNOGEN	Synthesized peptide derived from human Gab 2 around the non-phosphorylation site of S159.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat,Monkey
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.