

PR Monoclonal Antibody

OVERVIEW

CATALOG NO.	PTX17359	SYNONYMS
SIZE	100 ug	
HOST	Mouse	

PRODUCT INFORMATION

GENE NAME	PGR
PROTEIN NAME	Progesterone receptor
IMMUNOGEN	Synthetic Peptide

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	IHC IF
RECOMMENDED DILUTION	IHC: 1:200 IF: 1:100-200
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.