

Caspase-8 Monoclonal Antibody

OVERVIEW

CATALOG NO.	PTX17388	SYNONYMS
SIZE	100 ug	
HOST	Mouse	

PRODUCT INFORMATION

GENE NAME	CASP8
PROTEIN NAME	Caspase-8
IMMUNOGEN	Recombinant Protein

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB IHC
RECOMMENDED DILUTION	WB: 1:1000-2000 IHC:1:200-500
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.