

HSP90 α Monoclonal antibody

OVERVIEW

CATALOG NO.	PTX17490	SYNONYMS
SIZE	100 ug	HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86;
HOST	Rabbit	HSP86; Renal carcinoma antigen NY-REN-38

PRODUCT INFORMATION

GENE NAME	HSP90 α
PROTEIN NAME	HSP90 α
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human HSP90A. at AA range: 660-740

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.