

HSC70 Monoclonal antibody

OVERVIEW

CATALOG NO.	PTX17492	SYNONYMS
SIZE	100 ug	HSPA8; HSC70; HSP73; HSPA10; Heat shock cognate
HOST	Rabbit	71 kDa protein; Heat shock 70 kDa protein 8

PRODUCT INFORMATION

GENE NAME	HSC70
PROTEIN NAME	HSC70
IMMUNOGEN	Synthesized peptide derived from the Internal region of human HSC 70.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.