

HSP90 α Monoclonal antibody

OVERVIEW

| | | |
|--------------------|----------|--|
| CATALOG NO. | PTX17506 | SYNONYMS |
| SIZE | 100 ug | HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86; HSP86; Renal carcinoma antigen NY-REN-38 |
| HOST | Rabbit | |

PRODUCT INFORMATION

| | |
|---------------------|--|
| GENE NAME | HSP90 α |
| PROTEIN NAME | HSP90 α |
| IMMUNOGEN | Synthesized peptide derived from the C-terminal region of human HSP90A. at AA range: 660-740 |

APPLICATIONS AND STORAGE

| | |
|--------------------------------|--|
| SPECIES REACTIVITY | Human, Mouse, Rat |
| TESTED APPLICATIONS | WB |
| RECOMMENDED DILUTION | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications. |
| DELIVERY CONDITIONS | Blue Ice |
| STABILITY & STORAGE | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles |

For research use only.