

## MEK2 Mouse Monoclonal Antibody(2C3)

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX17510	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	MAP2K2; MEK2; MKK2; PRKMK2; Dual specificity
<b>HOST</b>	Rabbit	mitogen-activated protein kinase kinase 2; MAP kinase kinase 2; MAPKK 2; ERK activator kinase 2; MAPK/ERK kinase 2; MEK 2

### PRODUCT INFORMATION

---

#### GENE NAME

**PROTEIN NAME** Dual specificity mitogen-activated protein kinase kinase 2

**IMMUNOGEN** Synthesized peptide derived from human MEK-2 around the non-phosphorylation site of T394.

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human, Mouse, Rat
<b>TESTED APPLICATIONS</b>	WB
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. ELISA: 1/10000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**