

PI3 Kinase P85 α Mouse Monoclonal Antibody(2D2)

OVERVIEW

CATALOG NO.	PTX17514	SYNONYMS
SIZE	100 ug	PIK3R1; GRB1; Phosphatidylinositol 3-kinase regulatory subunit alpha; PI3-kinase regulatory subunit alpha; PI3K regulatory subunit alpha; PtdIns-3-kinase regulatory subunit alpha; Phosphatidylinositol 3-kinase 85 kDa regulatory subunit alpha; PI3-kinase s
HOST	Mouse	

PRODUCT INFORMATION

GENE NAME

PROTEIN NAME	PI3 Kinase P85 α
IMMUNOGEN	Recombinant Protein

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Rat, Mouse
TESTED APPLICATIONS	WB IHC
RECOMMENDED DILUTION	WB 1:1000-2000, IHC 1:100-200
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.