

Smad3 Mouse Monoclonal Antibody(8B5)

OVERVIEW

CATALOG NO.	PTX17516	SYNONYMS
SIZE	100 ug	SMAD3; MADH3; Mothers against decapentaplegic
HOST	Mouse	homolog 3; MAD homolog 3; Mad3; Mothers against DPP homolog 3; hMAD-3; JV15-2; SMAD family member 3; SMAD 3; Smad3; hSMAD3

PRODUCT INFORMATION

GENE NAME

PROTEIN NAME Smad3

IMMUNOGEN Recombinant Protein

APPLICATIONS AND STORAGE

SPECIES REACTIVITY Human, Mouse, Rat

TESTED APPLICATIONS WB IHC

RECOMMENDED DILUTION WB 1:1000-2000, IHC 1:100-200

DELIVERY CONDITIONS Blue Ice

STABILITY & STORAGE PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
4°C for short term (1 week), store at -20°C to -80°C for long term (1 year);
Avoid repeated freeze-thaw cycles

For research use only.