

STAT3 Mouse Monoclonal Antibody(3F4)

OVERVIEW

CATALOG NO.	PTX17517	SYNONYMS
SIZE	100 ug	STAT3; APRF; Signal transducer and activator of
HOST	Mouse	transcription 3; Acute-phase response factor

PRODUCT INFORMATION

GENE NAME

PROTEIN NAME STAT3

IMMUNOGEN Recombinant Protein

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB IHC
RECOMMENDED DILUTION	WB 1:1000-2000, IHC 1:100-200
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.