

GRP78/Bip Mouse Monoclonal Antibody(8G3)

OVERVIEW

CATALOG NO.	PTX17518	SYNONYMS
SIZE	100 ug	HSPA5; GRP78; 78 kDa glucose-regulated protein;
HOST	Mouse	GRP-78; Endoplasmic reticulum luminal Ca(2+)-binding protein grp78; Heat shock 70 kDa protein 5; Immunoglobulin heavy chain-binding protein; BiP

PRODUCT INFORMATION

GENE NAME

PROTEIN NAME GRP78/Bip

IMMUNOGEN Synthetic Peptide

APPLICATIONS AND STORAGE

SPECIES REACTIVITY Human, Rat

TESTED APPLICATIONS WB IHC

RECOMMENDED DILUTION WB 1:1000-2000, IHC 1:100-200

DELIVERY CONDITIONS Blue Ice

STABILITY & STORAGE PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
4°C for short term (1 week), store at -20°C to -80°C for long term (1 year);
Avoid repeated freeze-thaw cycles

For research use only.