

TBP/TATA Binding Protein Mouse Monoclonal Antibody(4H2)

OVERVIEW

CATALOG NO.	PTX17521	SYNONYMS
SIZE	100 ug	TBP; GTF2D1; TF2D; TFIID; TATA-box-binding protein; TATA sequence-binding protein; TATA-binding factor; TATA-box factor; Transcription initiation factor
HOST	Mouse	TFIID TBP subunit

PRODUCT INFORMATION

GENE NAME

PROTEIN NAME TBP/TATA Binding Protein

IMMUNOGEN Recombinant Protein

APPLICATIONS AND STORAGE

SPECIES REACTIVITY Human, Mouse, Rat

TESTED APPLICATIONS WB

RECOMMENDED DILUTION WB 1:1000-2000

DELIVERY CONDITIONS Blue Ice

STABILITY & STORAGE PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
4°C for short term (1 week), store at -20°C to -80°C for long term (1 year);
Avoid repeated freeze-thaw cycles

For research use only.