

Desmin Mouse Monoclonal Antibody(1B12)

OVERVIEW

CATALOG NO.	PTX17523	SYNONYMS
SIZE	100 ug	DES; Desmin
HOST	Mouse	

PRODUCT INFORMATION

GENE NAME

PROTEIN NAME	Desmin
IMMUNOGEN	Recombinant Protein

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	IHC
RECOMMENDED DILUTION	IHC 1:100-200
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.