

## Cyclophilin B Mouse Monoclonal Antibody(1B3)

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX17524	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	PPIB; CYPB; Peptidyl-prolyl cis-trans isomerase B;
<b>HOST</b>	Mouse	PPIase B; CYP-S1; Cyclophilin B; Rotamase B; S-cyclophilin; SCYLP

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	
<b>PROTEIN NAME</b>	Cyclophilin B
<b>IMMUNOGEN</b>	Synthetic Peptide

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human, Mouse, Rat
<b>TESTED APPLICATIONS</b>	WB IHC
<b>RECOMMENDED DILUTION</b>	WB 1:1000-2000, IHC 1:100-200
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**