

TBP/TATA Binding Protein Mouse Monoclonal Antibody(1F6)

OVERVIEW

CATALOG NO.	PTX17525	SYNONYMS
SIZE	100 ug	TBP; GTF2D1; TF2D; TFIID; TATA-box-binding protein; TATA sequence-binding protein; TATA-binding factor; TATA-box factor; Transcription initiation factor
HOST	Mouse	TFIID TBP subunit

PRODUCT INFORMATION

GENE NAME

PROTEIN NAME TBP/TATA Binding Protein

IMMUNOGEN Recombinant Protein

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	IHC
RECOMMENDED DILUTION	IHC 1:50-200
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.