

P44/42 MAPK(ERK1/2) Mouse Monoclonal Antibody(6B1)

OVERVIEW

CATALOG NO.	PTX17526	SYNONYMS
SIZE	100 ug	MAPK3; ERK1; PRKM3; Mitogen-activated protein kinase 3; MAP kinase 3; MAPK 3; ERT2; Extracellular signal-regulated kinase 1; ERK-1; Insulin-stimulated MAP2 kinase; MAP kinase isoform p44; p44-MAPK; Microtubule-associated protein 2 kinase; p
HOST	Mouse	

PRODUCT INFORMATION

GENE NAME	
PROTEIN NAME	P44/42 MAPK(ERK1/2)
IMMUNOGEN	Synthetic Peptide

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB IHC
RECOMMENDED DILUTION	WB 1:1000-2000, IHC 1:50-100
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.