

## **P44/42 MAPK(ERK1/2) Mouse Monoclonal Antibody(6C10)**

### **OVERVIEW**

---

<b>CATALOG NO.</b>	PTX17528	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	MAPK3; ERK1; PRKM3; Mitogen-activated protein kinase 3; MAP kinase 3; MAPK 3; ERT2; Extracellular signal-regulated kinase 1; ERK-1; Insulin-stimulated MAP2 kinase; MAP kinase isoform p44; p44-MAPK; Microtubule-associated protein 2 kinase; p
<b>HOST</b>	Mouse	

### **PRODUCT INFORMATION**

---

<b>GENE NAME</b>	
<b>PROTEIN NAME</b>	P44/42 MAPK(ERK1/2)
<b>IMMUNOGEN</b>	Synthetic Peptide

### **APPLICATIONS AND STORAGE**

---

<b>SPECIES REACTIVITY</b>	Human, Mouse, Rat
<b>TESTED APPLICATIONS</b>	WB IHC
<b>RECOMMENDED DILUTION</b>	WB 1:1000-2000, IHC 1:50-100
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**