

α -tubulin (Acetyl Lys40) Mouse Monoclonal Antibody(4A8)

OVERVIEW

CATALOG NO.	PTX17529	SYNONYMS
SIZE	100 ug	
HOST	Mouse	

PRODUCT INFORMATION

GENE NAME

PROTEIN NAME α -tubulin (Acetyl Lys40)

IMMUNOGEN Synthetic Peptide

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB IHC
RECOMMENDED DILUTION	WB 1:1000-2000, IHC 1:50-100
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.