

I κ B β Mouse Monoclonal Antibody(1F3)

OVERVIEW

CATALOG NO.	PTX17531	SYNONYMS
SIZE	100 ug	NFKBIB; IKBB; TRIP9; NF-kappa-B inhibitor beta; NF-kappa-BIB; I-kappa-B-beta; I κ B-B; I κ B-beta;
HOST	Mouse	IkappaBbeta; Thyroid receptor-interacting protein 9; TR-interacting protein 9; TRIP-9

PRODUCT INFORMATION

GENE NAME

PROTEIN NAME I κ B β

IMMUNOGEN Recombinant Protein

APPLICATIONS AND STORAGE

SPECIES REACTIVITY Human, Mouse, Rat

TESTED APPLICATIONS IHC

RECOMMENDED DILUTION IHC 1:100-200

DELIVERY CONDITIONS Blue Ice

STABILITY & STORAGE PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
4°C for short term (1 week), store at -20°C to -80°C for long term (1 year);
Avoid repeated freeze-thaw cycles

For research use only.