

HSPB8/HSP22 Mouse Monoclonal Antibody(2C3)

OVERVIEW

CATALOG NO.	PTX17535	SYNONYMS
SIZE	100 ug	HSPB8; CRYAC; E2IG1; HSP22; PP1629; Heat shock protein beta-8; HspB8; Alpha-crystallin C chain; E2-induced gene 1 protein; Protein kinase H11; Small stress protein-like protein HSP22
HOST	Mouse	

PRODUCT INFORMATION

GENE NAME

PROTEIN NAME HSPB8/HSP22

IMMUNOGEN Recombinant Protein

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB
RECOMMENDED DILUTION	WB 1:1000-2000
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.