

MAO-A Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17644	SYNONYMS
SIZE	100 ug	MAOA; Amine oxidase [flavin-containing] A;
HOST	Rabbit	Monoamine oxidase type A; MAO-A

PRODUCT INFORMATION

GENE NAME	MAOA
PROTEIN NAME	Amine oxidase [flavin-containing] A
IMMUNOGEN	Synthesized peptide derived from the Internal region of human MAO-A.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.