

IL-1F5 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17658	SYNONYMS
SIZE	100 ug	IL36RN; FIL1D; IL1F5; IL1HY1; IL1L1; IL1RP3;
HOST	Rabbit	Interleukin-36 receptor antagonist protein; FIL1 delta; IL-1-related protein 3; IL-1RP3; Interleukin-1 HY1; IL- 1HY1; Interleukin-1 delta; IL-1 delta; Interleukin-1 family member 5; IL-1F5; Interleukin-1 recepto

PRODUCT INFORMATION

GENE NAME	IL36RN
PROTEIN NAME	Interleukin-36 receptor antagonist protein
IMMUNOGEN	Synthesized peptide derived from the Internal region of human IL-1F5.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.