

FOG-2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17669	SYNONYMS
SIZE	100 ug	ZFPM2; FOG2; ZNF89B; Zinc finger protein ZFPM2;
HOST	Rabbit	Friend of GATA protein 2; FOG-2; Friend of GATA 2; hFOG-2; Zinc finger protein 89B; Zinc finger protein multitype 2

PRODUCT INFORMATION

GENE NAME	ZFPM2
PROTEIN NAME	Zinc finger protein ZFPM2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human FOG-2.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse
TESTED APPLICATIONS	WB IHC ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.