

Synuclein- α Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17684	SYNONYMS
SIZE	100 ug	SNCA; NACP; PARK1; Alpha-synuclein; Non-A beta
HOST	Rabbit	component of AD amyloid; Non-A4 component of amyloid precursor; NACP

PRODUCT INFORMATION

GENE NAME	SNCA
PROTEIN NAME	Alpha-synuclein
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Synuclein- α .

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB IHC ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.