

## **FITC\* Goat Anti-Mouse IgG (H+L)**

### **OVERVIEW**

---

<b>CATALOG NO.</b>	PTX17818-300
<b>SIZE</b>	300ul
<b>HOST</b>	Goat

### **PRODUCT INFORMATION**

---

<b>CONCENTRATION</b>	~4mg/mL
<b>CLONALITY</b>	Polyclonal Antibody

### **APPLICATIONS & STORAGE**

---

<b>REACTIVITY</b>	Mouse
<b>APPLICATIONS</b>	WB IHC IP IF FCM ELISA
<b>DILUTION</b>	WB 1:5000-20000, IHC 1:200-2000, IF 1:200-2000, IP 1:200-1:2000, IP 1:200-2000, FCM 1:200-2000, ELISA 1:5000-20000
<b>FORMULATION</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>PURIFICATION</b>	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
<b>DELIVERY CONDITION</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles