

BIOTIN* Goat Anti-Rabbit IgG (H+L)

OVERVIEW

CATALOG NO.	PTX17821-300
SIZE	300ul
HOST	Goat

PRODUCT INFORMATION

CONCENTRATION	~4mg/mL
CLONALITY	Polyclonal Antibody

APPLICATIONS & STORAGE

REACTIVITY	Rabbit
APPLICATIONS	WB IHC IP IF FCM ELISA
DILUTION	WB 1:5000-20000, IHC 1:200-2000, IF 1:200-2000, IP 1:200-1:2000, IP 1:200-2000, FCM 1:200-2000, ELISA 1:5000-20000
FORMULATION	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
PURIFICATION	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
DELIVERY CONDITION	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles