

Di-Methyl-Histone H3 (K27) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17824	SYNONYMS
SIZE	100 ug	
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	HIST1H3A/HIST1H3/HIST1H3C/HIST1H3D/HIST1H3E/HIST1H3F/HIST1H3G/HIST1H3H/HIST1H3I/HIST1H3J/HIST2H3A/HIST2H3C/HIST2H3D/H3F3A/H3F3B/H3F3C
PROTEIN NAME	Histone H3.1/Histone H3.2/Histone H3.3/Histone H3.3C
IMMUNOGEN	Synthesized peptide derived from human Histone H3 around the methylation site of K27.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.