

Acetyl-14-3-3-pan (K51/49) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17826	SYNONYMS
SIZE	100 ug	
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	YWHAB/YWHAG/YWHAQ/YWHAZ/SFN
PROTEIN NAME	14-3-3 protein beta/alpha/14-3-3 protein gamma/14-3-3 protein theta/14-3-3 protein zeta/delta/14-3-3 protein sigma
IMMUNOGEN	Synthesized peptide derived from human 14-3-3-pan around the acetylation site of K51.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.