

Acetyl-p300 (K1542) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17827	SYNONYMS
SIZE	100 ug	
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	EP300
PROTEIN NAME	Histone acetyltransferase p300
IMMUNOGEN	Synthesized peptide derived from human p300 around the acetylation site of K1542.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.