

Acetyl-Stat3 (K87) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX17830	SYNONYMS
SIZE	100 ug	
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	STAT3
PROTEIN NAME	Signal transducer and activator of transcription 3
IMMUNOGEN	Synthesized peptide derived from human Stat3 around the acetylation site of K87.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.