

## Mouse Monoclonal antibody to Rat EP24.15

### OVERVIEW

---

<b>CATALOG NO.</b>	PX-A1050-1H6-2G2-L
<b>SIZE</b>	100µL
<b>CLONALITY</b>	Monoclonal Antibody
<b>HOST</b>	Mouse
<b>CLONE NAME</b>	1H6-2G2
<b>ISOTYPE</b>	IgG1, y

### APPLICATIONS AND STORAGE

---

<b>REACTIVITY</b>	Human, Mouse, Rat, Primate, Rabbit
<b>APPLICATIONS</b>	WB, ICC/IF, IHC, IHC-P
<b>FORM</b>	liquid/frozen
<b>DELIVERY CONDITION</b>	Dry Ice
<b>STABILITY &amp; STORAGE</b>	Store at -20°C. Aliquot to avoid repeated freezing and thawing.