

Mouse Monoclonal antibody to Rat EP24.15

OVERVIEW

CATALOG NO.	PX-A1050-4C2-2B10-L
SIZE	100µL
CLONALITY	Monoclonal Antibody
HOST	Mouse
CLONE NAME	4C2-2B10
ISOTYPE	IgG2a, k

APPLICATIONS AND STORAGE

REACTIVITY	Human, Mouse, Rat, Primate, Rabbit
APPLICATIONS	WB, ICC/IF, IHC, IHC-P
FORM	liquid/frozen
DELIVERY CONDITION	Dry Ice
STABILITY & STORAGE	Store at -20°C. Aliquot to avoid repeated freezing and thawing.