

Mouse Monoclonal antibody to Rat EP24.15

OVERVIEW

| | |
|--------------------|---------------------|
| CATALOG NO. | PX-A1050-4C2-2B9-L |
| SIZE | 100µL |
| CLONALITY | Monoclonal Antibody |
| HOST | Mouse |
| CLONE NAME | 4C2-2B9 |
| ISOTYPE | IgG2a, k |

APPLICATIONS AND STORAGE

| | |
|--------------------------------|---|
| REACTIVITY | Human, Mouse, Rat, Primate, Rabbit |
| APPLICATIONS | WB, ICC/IF, IHC, IHC-P |
| FORM | liquid/frozen |
| DELIVERY CONDITION | Dry Ice |
| STABILITY & STORAGE | Store at -20°C. Aliquot to avoid repeated freezing and thawing. |