

MRP63-mitochondrial ribosomal protein 63 Gene

Reference: PTXBC000002

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MRP63; bMRP63; ribosomal protein 63, mitochondrial; hMRP63; mitochondrial ribosomal protein 63; mitochondrial ribosomal protein bMRP63; mitochondrial ribosomal protein L57		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein 63
GENE ID	78988
ORIGIN SPECIES	Human
ACCESSIONS	BC000002
VECTOR	pDONR233
NCBI SYMBOL	MRP63
LENGTH	309 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgttcctgactgcgctcctctggcggcggccgattcccgccgtcagtgatcgggaagcaccggcggcgggttcgtgctgtgcgccaagcagaacatgatccgccgctggagatgagcgggagaaccattactggctgagcatgcc  
ctacatgacccgggagcaggagcggcggccacggcgggtgcaggaggaggccttcgaggccataaaggcggccgccactccaagttccccccgatagattcattgcggaccagctcgacctctcaatgtccaagaatgtcctaa
```