

GTPBP8-GTP-binding protein 8 (putative) Gene

Reference: PTXBC000003

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC135; GTP-binding protein 8; GTP binding protein 8 (putative)		

For research use only.

GENE DESCRIPTION	GTP-binding protein 8 (putative)
GENE ID	29083
ORIGIN SPECIES	Human
ACCESSIONS	BC000003
VECTOR	pDONR233
NCBI SYMBOL	GTPBP8
LENGTH	513 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggcggcggcgggctgcggctgggagcgggaaagactcttfgaaatgcctcgggtctagagc gactgagccgctataatgcacgtcccaagcttttctgaggtgctgctgctgccgaagcagcagctgaggaagctgctgtaccgctgca
ggaaagtagagcgggttcctcgcctccctacgggagcgaagacctcacctgcgtatcttfgaccnaagcccggaggacatagccaggcgggacaacatcttcacggccactgaaccgaaccgatcgcactacgtcagctccgctccgtatcgacca
cgccccgacctcccgccagaggtgttttataggcagaagcaatgttgaaaatcatctctaatcaagccttttttctactgccccctgaggtgaaagtcagagctcctcaaaaaccagacacacaaaataatgtttttcaagtggaaa
acatfttacagctggtggacatgccaggttatgctttag
```