

PRELID1-PRELI domain containing 1 Gene

Reference: PTXBC000007

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-106; PX19; SBBI12; PRELI domain-containing protein 1, mitochondrial; 25 kDa protein of relevant evolutionary and lymphoid interest; px19-like protein; PRELI domain containing 1		

For research use only.

GENE DESCRIPTION	PRELI domain containing 1
GENE ID	27166
ORIGIN SPECIES	Human
ACCESSIONS	BC000007
VECTOR	pDONR233
NCBI SYMBOL	PRELID1
LENGTH	660 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggtgaagtattctctggccagagcgtgctcggagttcctgggaccaagtgttcgccctcttggcagcgggtaccggaatccctatagcaaacatgtcttgacggaagacatagtacaccgggaggtgaccctgaccagaaactgctgccccg
gcgactcctgaccaagaccaacaggatgccacgctgggccgagcgaactttcctgccaatgttctcactcgggtgactcctggaggactctattgtggaccacagaatcagaccatgactaccttcaactggaacatcaaccacgcccggctgat
gggtggagggaacgatgtgttactgtggaactctgacaacagtggctggactgaaatccgccgggaagcctgggtctcctctagcttattgtgtctccagagctgccaggaattgtcttcccgggttcaaaaagcaacgtgaccaagactatga
aggggtttgaatatactctgctaagctgcaaggcgaaggcccttccaaaacactgtgtgagacagccaaggaaagcaaggagaaaggcaaggacggcactggcagctacagagaaggccaaggacctgccagcaaggcggccaccaaga
agcagcagcagcagcaacagttttag
```