

NHP2-NHP2 ribonucleoprotein homolog (yeast) Gene

Reference: PTXBC000009

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NHP2 ribonucleoprotein; snoRNP protein NHP2; NHP2-like protein; NHP2 ribonucleoprotein homolog; DKCB2; NHP2P; NOLA2; H/ACA ribonucleoprotein complex subunit 2; nucleolar protein family A, member 2 (H/ACA small nucleolar RNPs)		

For research use only.

GENE DESCRIPTION	NHP2 ribonucleoprotein homolog (yeast)
GENE ID	55651
ORIGIN SPECIES	Human
ACCESSIONS	BC000009
VECTOR	pDONR223
NCBI SYMBOL	NHP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgacaaaataaaggcagatcccgcacggcccgaggctcagcggagcgtgttcggggagcgcacctaccaggagctgtggtcaaccagaaccccatcgcgcagccctggtctcgcgcctcacgcggaagctctacaatgcatcaagaagcgggtgaagcagaagc
agattcggcgcgggtgaaagaggttcagaaattgtcaacaaggagaaaaaggatcatgttttggcaggagacacacagccattgaggtatactgccatctcccagtcattgtgtgaggaccgaaatttgcctatgtctatccccttaagacggacctgggtgcagccgcaggc
tccaagcggccacctgtgtgataatggtcaagcccatgaggagtaccaggagcttacgatgagtgctggaggaggtgcagtcctgcccctaccctatga
```