

ARPC5L-actin related protein 2/3 complex, subunit 5-like Gene

Reference: PTXBC000018

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARC16-2; actin-related protein 2/3 complex subunit 5-like protein; arp2/3 complex 16 kDa subunit 2; actin related protein 2/3 complex subunit 5 like		

For research use only.

GENE DESCRIPTION	actin related protein 2/3 complex, subunit 5-like
GENE ID	81873
ORIGIN SPECIES	Human
ACCESSIONS	BC000018
VECTOR	pDONR223
NCBI SYMBOL	ARPC5L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccgaacacgctgtcctcgcctccgccgggtggacatcgacaatttgacgagacaatttggacgagcagaggcggggcggcgggcggcgagccaggccccgaccgagcaggtggacggcctcgcggcaaggggacatgcttcggccattccat
gcagcctgcgggaactcctccgtcaacccaagaatcaagctgtgaaggagcagcccaggcgtggtgctgaaagtgtctcacaactcgaagcagtgcagtgcagcaggcgtgcagtgcagtcactggacagaaacggcgtgactgtaatgaagtacattataaaggcttgagaagc
ccacagaaaatagcagcagtgcttactcagtggcacgaaaaggccttagcagtaggagtaggactaggtccattataagagtcttacagcaagaagactgttt