

RPS19-ribosomal protein S19 Gene

Reference: PTXBC000023

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	DBA; DBA1; S19; 40S ribosomal protein S19; ribosomal protein S19		

For research use only.

GENE DESCRIPTION	ribosomal protein S19
GENE ID	6223
ORIGIN SPECIES	Human
ACCESSIONS	BC000023
VECTOR	pDONR223
NCBI SYMBOL	RPS19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctggagtactgtaaaagacgtgaaccagcaggagttcgtcagagctctggcagcctctcaaaaagtcgggaagctgaaagtccccgaatgggtgalaccgtcaagctggccaagcacaagaagctgtcctctacgatgagaactggtctacacgcgagctgctccaca
gcgcgacactgtactccgggtggcgtgggtgtgctccatgaccaagatctatggggacgtcagagaaaacggctcatgccagccactcagccgaggtccaagagtggtgccgccgggtctccaagccctggagggctgaaaatggtgaaaaggaccaagatg
gcggccgcaaacgtgacacctcagggacaagagatctggacagaatgccggacaggtggcagctgccaacaagaagcattag