

FAM60A-family with sequence similarity 60, member A Gene

Reference: PTXBC000024

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM60A; C12orf14; TERA; tera protein homolog; family with sequence similarity 60 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 60, member A
GENE ID	58516
ORIGIN SPECIES	Human
ACCESSIONS	BC000024
VECTOR	pDONR233
NCBI SYMBOL	FAM60A
LENGTH	666 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgtttggtttcacaagccaagatgtaccgaagtatagagggctgctgtattgcagagctaaagtctccagttctgattcactgacagtaaacgctatgaaaaggactccagagctgttttgattgcatgagactcgtcaggagacatctgcaatgc
ctgtgctcgtgctgtgaaaagatggaagaagtgtccagcagggatcaaaaaaaaaaactggaatcatgtggtatgcaagggctggaccagctctaaagactacattgaaaccaaaagaaagtgaaactctatctgggaacaggataaaaaagcaaccag
atcagtaaacctcagaaggaaltnaacgtcataattctgatgctcagctaccacctcaagtgcctccccagctcaatctcctgtttacagtaaccagtcagatgacggctcagatacagagatggcttctgttctaacagaacaccagtttttctttta
gatctcactactggaagacagagaagatgttgggatcatataaaggccgtttggggaagtcctcattgacacacatctctcaagccttctgctgagcaataagaagcagctgctgagaagccagaggagcaggggcccagagcctctgcc
atctccactcaggagtgtga
```