

TMED3-transmembrane emp24 protein transport domain containing 3 Gene

Reference: PTXBC000027

| | | | |
|----------------|---|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | C15orf22; P24B; p24g4; p26; transmembrane emp24 domain-containing protein 3; integral type I protein; membrane protein p24B; p24 family protein gamma-4; p24gamma4; testis tissue sperm-binding protein Li 48e; transmembrane emp24 protein transport domain co | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | transmembrane emp24 protein transport domain containing 3 |
| GENE ID | 23423 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC000027 |
| VECTOR | pDONR233 |
| NCBI SYMBOL | TMED3 |
| LENGTH | 654 bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

```
atgggcagcactgtcccgcctccctccctcgtctgcttctgctgctgctcctgcccgggcccagcagccctcggggccgagctcacttccgagctgcccgacaacccaagcagctcctccacgaggagggtgagcagggcgtgaagttctc
cctggattaccagggtcactcactggaggcactacgatggtgactgctatgtagaggacccccaggggaacaccatctacagagaacgaagaagcagtagcacagcttcacgtaccgggctgaagtaagggcgtttatcagtttgctcagtaatga
gtttccaccittctcacaagaccgtctacttgaacttcaagtggcgatgagcctccattctcccagacatggggaacagggtcacagctctcaccagatggagtccctcgtgacctccatgaggctctgaaaacgggtgaltgactcccaga
cgcattaccggctcgggagcccaggaccgggcccagcgggaagacctaatagccgagctcttactgtctgttggcgagacgattgccctgtctgtggtcagctcagtcaggtgctactgttgaagaactcttcacagaaaacgacctca
gcaggcagtcactcactag
```