

TMED3-transmembrane emp24 protein transport domain containing 3 Gene

Reference: PTXBC000027

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C15orf22; P24B; p24g4; p26; transmembrane emp24 domain-containing protein 3; integral type I protein; membrane protein p24B; p24 family protein gamma-4; p24gamma4; testis tissue sperm-binding protein Li 48e; transmembrane emp24 protein transport domain co		

For research use only.

GENE DESCRIPTION	transmembrane emp24 protein transport domain containing 3
GENE ID	23423
ORIGIN SPECIES	Human
ACCESSIONS	BC000027
VECTOR	pDONR223
NCBI SYMBOL	TMED3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgggcagcactgtcccgcctccgctcctctgtctctgctcctctgcccggccgagcagcctctgcccggccgagctcacttcgagctgccggacaacccaagcagtgctccacagagaggtggagcaggcgctgaagtctccctggattaccaggtcactc  
gagccactacgatgttactgctatgagaggaccocagggaacaccatctacagagaaacgaagaagcagtacgacgctcactcaccggctgaagtcaagggctttatcagtttgcctcagtaatgagtttccacctctctcaaacaccgctactttgacttcaagggg  
cgatgagcctccattctcccagacatggggaacagggtcacagctctaccagatggagtcgctgacatccatgaggtctgaaaacggtgattgactcccagacgattaccggctcggggaggcccaggaccggcccagcgggaagacitaaagccgagctct  
tactgtctgtggcagacattgccctgtctgttcagctcagtcaggtgctactgttgaagaacctctcacagaaaaacacccatcagcaggcagtcactcctag
```